

S.C. CLINI-LAB S.R.L

RO-540342, Tîrgu-Mureş, str. Rodnei nr. 15, tel: 0265-218.784, fax:0265-218.407

RO-010917, Bucureşti, str. C-tin Bălăcescu nr. 5, tel: 021 - 313.93.55, fax: 021 - 313.93.54

web: www.clinilab.ro, email: office@clinilab.ro

DISCURI PENTRU ANTIBIOGRAMĂ BIOLAB

COD	ANTIBIOTIC	COD2	ANTIBIOTIC2
Conţinut: cartuş cu 50 discuri			
AK-10	Amikacin, 10 µg	CTC-40	Cefotaxime + Clavulanic acid, 30+10 µg
AK-30	Amikacină, 30 µg	CTX-5	Cefotaxime, 5 µg
AN-60	Aminosidin, 60 µg	CTX-10	Cefotaxime, 10 µg
AMC-3	Amoxicillin + Clavulanic acid, 2+1 µg	CTX-30	Cefotaxime, 30 µg
AMC-30	Amoxicillin + Clavulanic acid, 20+10 µg	CTT-30	Cefotetan, 30 µg
AX-2	Amoxicillin, 2 µg	CTM-30	Cefotiam, 30 µg
AX-10	Amoxicillin, 10 µg	FOX-10	Cefoxitin, 10 µg
AX-25	Amoxicillin, 25 µg	FOX-30	Cefoxitin, 30 µg
AX-30	Amoxicillin, 30 µg	CPO-30	Cefpirome, 30 µg
AMB-100	Amphotericin B, 100 U	CPC-11	Cefpodoxime + Clavulanic acid, 10+1 µg
APX-30	Ampicillin + Cloxacillin, 25+5 µg	CPD-10	Cefpodoxime, 10 µg
SAM-20	Ampicillin + Sulbactam, 10+10 µg	CPR-30	Cefprozil, 30 µg
SAM-30	Ampicillin + Sulbactam, 20+10 µg	CFS-30	Cefsulodin, 30 µg
AM-2	Ampicillin, 2 µg	CZC-40	Ceftazidime + Clavulanic acid, 30+10 µg
AM-10	Ampicillin, 10 µg	CAZ-10	Ceftazidime, 10 µg
AM-25	Ampicillin, 25 µg	CAZ-30	Ceftazidime, 30 µg
AM-30	Ampicillin, 30 µg	CBT-30	Ceftibuten, 30 µg
APR-15	Apramycin, 15 µg	FUR-30	Ceftiofur, 30 µg
APR-30	Apramycin, 30 µg	ZOX-30	Ceftizoxime, 30 µg
AZM-15	Azithromycin, 15 µg	CRO-5	Ceftriaxone, 5 µg
AZM-30	Azithromycin, 30 µg	CRO-10	Ceftriaxone, 10 µg
AZL-75	Azlocillin, 75 µg	CRO-30	Ceftriaxone, 30 µg
ATM-30	Aztreonam, 30 µg	CXM-5	Cefuroxime, 5 µg
B-10	Bacitracin, 10 U	CXM-30	Cefuroxime, 30 µg
APB-250	Boronic acid, 250 µg	CXA-30	Cefuroxime axetil, 30 µg
APB-300	Boronic acid, 300 µg	CL-30	Cephalexin, 30 µg
APB-400	Boronic acid, 400 µg	KF-30	Cephalothin, 30 µg
PY-25	Carbenicillin, 25 µg	CE-30	Cephradine, 30 µg
PY-100	Carbenicillin, 100 µg	C-5	Chloramphenicol, 5 µg
CEC-30	Cefaclor, 30 µg	C-10	Chloramphenicol, 10 µg
CFR-30	Cefadroxil, 30 µg	C-30	Chloramphenicol, 30 µg
MA-30	Cefamandole, 30 µg	C-50	Chloramphenicol, 50 µg

CZ-30	Cefazolin, 30 µg	CHT-30	Chlortetracycline, 30 µg
CD-5	Cefdinir, 5 µg	CIN-100	Cinoxacin, 100 µg
CDN-5	Cefditoren, 5 µg	CIP-1	Ciprofloxacin, 1 µg
FEP-30	Cefepime, 30 µg	CIP-5	Ciprofloxacin, 5 µg
CFM-5	Cefixime, 5 µg	CIP-10	Ciprofloxacin, 10 µg
CDZ-30	Cefodizime, 30 µg	CIP-30	Ciprofloxacin, 30 µg
CES-105	Cefoperazone + Sulbactam, 75+30 µg	CLR-2	Clarithromycin, 2 µg
CEP-30	Cefoperazone, 30 µg	CLR-5	Clarithromycin, 5 µg
CEP-75	Cefoperazone, 75 µg	CLR-15	Clarithromycin, 15 µg
DA-2	Clindamycin, 2 µg	CN-200	Gentamicin, 200 µg
DA-10	Clindamycin, 10 µg	GRS-10	Griseofulvin, 10 µg
CLT-10	Clotrimazole, 10 µg	IPM-10	Imipenem, 10 µg
CLT-300	Clotrimazole, 300 µg	IMC-20	Imipenem + Cilastatin, 10+10 µg
CX-1	Cloxacillin, 1 µg	ISP-30	Isepamicin, 30 µg
CX-5	Cloxacillin, 5 µg	INH-1	Isoniazid, 1 µg
CX-10	Cloxacillin, 10 µg	INH-5	Isoniazid, 5 µg
CX-30	Cloxacillin, 30 µg	ITR-10	Itraconazole, 10 µg
CT-10	Colistin, 10 µg	ITR-50	Itraconazole, 50 µg
CT-25	Colistin, 25 µg	K-5	Kanamycin, 5 µg
CT-50	Colistin, 50 µg	K-30	Kanamycin, 30 µg
CT-100	Colistin, 100 µg	K-1000	Kanamycin, 1000 µg
DAP-30	Daptomycin, 30 µg	KTC-10	Ketoconazole, 10 µg
DB-10	Dibekacin, 10 µg	KTC-50	Ketoconazole, 50 µg
DO-10	Doxycycline, 10 µg	LAT-30	Latamoxef, 30 µg
DO-30	Doxycycline, 30 µg	LEV-1	Levofloxacin, 1 µg
ECO-10	Econazole, 10 µg	LEV-5	Levofloxacin, 5 µg
EN-10	Enoxacin, 10 µg	LCN-75	Lincomycin + Neomycin, 15+60 µg
ENR-5	Enrofloxacin, 5 µg	LCS-109	Lincomycin + Spectinomycin, 9+100 µg
ENR-10	Enrofloxacin, 10 µg	L-2	Lincomycin, 2 µg
ETP-10	Ertapenem, 10 µg	L-10	Lincomycin, 10 µg
E-5	Erythromycin, 5 µg	L-15	Lincomycin, 15 µg
E-10	Erythromycin, 10 µg	LNZ-10	Linezolid, 10 µg
E-15	Erythromycin, 15 µg	LNZ-30	Linezolid, 30 µg
E-30	Erythromycin, 30 µg	LOM-10	Lomefloxacin, 10 µg
FFC-30	Florfenicol, 30 µg	LOM-30	Lomefloxacin, 30 µg
FL-5	Flucloxacillin, 5 µg	LOR-30	Loracarbef, 30 µg
FLU-10	Fluconazole, 10 µg	MEC-10	Mecillinam, 10 µg
FLU-25	Fluconazole, 25 µg	MEC-25	Mecillinam, 25 µg
AFY-1	Flucytosine, 1 µg	MEM-10	Meropenem, 10 µg
FLM-30	Flumequine, 30 µg	ME-5	Methicillin, 5 µg
FF-50	Fosfomicin, 50 µg	ME-10	Methicillin, 10 µg
FF-200	Fosfomicin, 200 µg	ME-30	Methicillin, 30 µg
FM-100	Framycetin, 100 µg	MET-2.5	Metronidazole, 2.5 µg
FX-15	Furazolidone, 15 µg	MET-4	Metronidazole, 4 µg
FX-50	Furazolidone, 50 µg	MET-5	Metronidazole, 5 µg
FX-100	Furazolidone, 100 µg	MET-50	Metronidazole, 50 µg
FA-5	Fusidic acid, 5 µg	MET-80	Metronidazole, 80 µg
FA-10	Fusidic acid, 10 µg	MEZ-30	Mezlocillin, 30 µg

FA-30	Fusidic acid, 30 µg	MEZ-75	Mezlocillin, 75 µg
FA-50	Fusidic acid, 50 µg	MCZ-10	Miconazole, 10 µg
GAT-5	Gatifloxacin, 5 µg	MI-25	Minocycline, 25 µg
GAT-10	Gatifloxacin, 10 µg	MI-30	Minocycline, 30 µg
GAT-30	Gatifloxacin, 30 µg	MOX-30	Moxalactam, 30 µg
CN-10	Gentamicin, 10 µg	MXF-1	Moxifloxacin, 1 µg
CN-30	Gentamicin, 30 µg	MXF-5	Moxifloxacin, 5 µg
CN-50	Gentamicin, 50 µg	MUP-5	Mupirocin, 5 µg
CN-120	Gentamicin 120 µg	MUP-20	Mupirocin, 20 µg
MUP-200	Mupirocin, 200 µg	PT-15	Pristinamycin, 15 µg
NAF-1	Nafcillin, 1 µg	SYN-15	Quinupristin + Dalfopristin, 4.5+10.5 µg
NA-30	Nalidixic acid, 30 µg	RA-2	Rifampin, 2 µg
N-5	Neomycin, 5 µg	RA-5	Rifampin, 5 µg
N-10	Neomycin, 10 µg	RA-30	Rifampin, 30 µg
N-30	Neomycin, 30 µg	RF-30	Rifamycin SV, 30 µg
NET-10	Netilmicin, 10 µg	RFX-40	Rifaximin, 40 µg
NET-30	Netilmicin, 30 µg	ROX-30	Roxithromycin, 30 µg
F-50	Nitrofurantoin, 50 µg	SIS-10	Sisomicin, 10 µg
F-100	Nitrofurantoin, 100 µg	SPS	Sodium Polyanethol Sulfonate
F-200	Nitrofurantoin, 200 µg	SPX-5	Sparfloxacin, 5 µg
F-300	Nitrofurantoin, 300 µg	SPX-10	Sparfloxacin, 10 µg
NR-100	Nitrofurazone, 100 µg	SPT-10	Spectinomycin, 10 µg
NTX-30	Nitroxoline, 30 µg	SPT-25	Spectinomycin, 25 µg
NOR-2	Norfloxacin, 2 µg	SPT-100	Spectinomycin, 100 µg
NOR-5	Norfloxacin, 5 µg	SP-100	Spiramycin, 100 µg
NOR-10	Norfloxacin, 10 µg	SP-30	Spiramycin, 30 µg
NV-30	Novobiocin, 30 µg	S-10	Streptomycin, 10 µg
NY-100	Nystatin, 100 U	S-25	Streptomycin, 25 µg
OFX-2	Ofloxacin, 2 µg	S-300	Streptomycin, 300 µg
OFX-5	Ofloxacin, 5 µg	SD-100	Sulfadiazine, 100 µg
OFX-10	Ofloxacin, 10 µg	SD-300	Sulfadiazine, 300 µg
OL-15	Oleandomycin, 15 µg	STZ-300	Sulfathiazole, 300 µg
OX-1	Oxacillin, 1 µg	ST-30	Sulfisoxazole, 30 µg
OX-5	Oxacillin, 5 µg	ST-300	Sulfisoxazole, 300 µg
OX-10	Oxacillin, 10 µg	SMZ-25	Sulphamethoxazole, 25 µg
OA-2	Oxolinic acid, 2 µg	SMZ-50	Sulphamethoxazole, 50 µg
T-30	Oxytetracycline, 30 µg	SMZ-100	Sulphamethoxazole, 100 µg
PEF-5	Pefloxacin, 5 µg	TEC-30	Teicoplanin, 30 µg
P-1	Penicillin, 1 U	TEL-15	Telithromycin, 15 µg
P-1.5	Penicillin, 1.5 U	TE-10	Tetracycline, 10 µg
P-2	Penicillin, 2 U	TE-25	Tetracycline, 25 µg
P-5	Penicillin, 5 U	TE-30	Tetracycline, 30 µg
P-10	Penicillin, 10 U	TP-30	Thiamphenicol, 30 µg
PI-20	Pipemidic acid, 20 µg	TIM-85	Ticarillin + Clavulanic acid, 75+10 µg
PI-30	Pipemidic acid, 30 µg	TIC-75	Ticarillin, 75 µg
TPZ-36	Piperacillin + Tazobactam, 30+6 µg	TGC-15	Tigecycline, 15 µg
TPZ-40	Piperacillin + Tazobactam, 30+10 µg	TIL-15	Tilmicosin, 15 µg

TPZ-85	Piperacillin + Tazobactam, 75+10 µg	TOB-10	Tobramycin, 10 µg
TPZ-110	Piperacillin + Tazobactam, 100+10 µg	TOB-30	Tobramycin, 30 µg
PRL-30	Piperacillin, 30 µg	TS-30	Trimethoprim + Sulfadiazine, 5+25 µg
PRL-75	Piperacillin, 75 µg	SXT-25	Trimethoprim+Sulphamethoxazole, 1.25+23.75 µg
PRL-100	Piperacillin, 100 µg	TMP-1.25	Trimethoprim, 1.25 µg
PIR-2	Pirlimycin, 2 µg	TMP-2.5	Trimethoprim, 2.5 µg
PB-50	Polymixin B, 50 U	TMP-5	Trimethoprim, 5 µg
PB-100	Polymixin B, 100 U	TMP-10	Trimethoprim, 10 µg
PB-300	Polymixin B, 300 U	TMP-25	Trimethoprim, 25 µg
TMP-30	Trimethoprim, 30 µg	VA-5	Vancomycin, 5 µg
S3-300	Triple Sulphonamides, 300 µg	VA-10	Vancomycin, 10 µg
TY-15	Tylosin, 15 µg	VA-30	Vancomycin, 30 µg
TY-30	Tylosin, 30 µg	VOR-1	Voriconazole, 1 µg
TY-60	Tylosin, 60 µg	X	X Factor Disc
V	V Factor Disc	X+V	X+V Factor Disc

Discuri pentru sensibilitate, diferențiale

B04	Bacitracină, 0,04 U	NV	Novobiocin, 5 µg
OP	Optochină, 5 µg		

Interpretarea rezultatelor este conform NCCLS